


SCULPO RETRAKTOREN MIT SILIKONHANDGRIFF
SCULPO RETRACTORS WITH SILIKON HANDLE




Der einzigartige Handgriff aus sterilisierbarem Silikon ermöglicht ein absolut ermüdungsfreies Arbeiten. Die anatomische Handform garantiert einen guten Halt und vermeidet ein Abrutschen selbst unter schwersten Bedingungen. Der Ø 5 mm Kaltlichtkanal mit Universaladapter hat eine äußerst hohe Lichtleiter-Qualität und leuchtet das komplette Operationsfeld bestmöglich aus. Der Kaltlichtkanal wird durch den Handgriff geführt wodurch ein OP-Feld entsteht und keine Kabel die Sicht beeinträchtigen.

The uniquely designed handle made of sterilizable silicon facilitates working without signs of fatigue. The anatomical design handle ensures a secure grip and avoids slippage under hard working conditions. The Ø 5 mm fiber optic channel with an universal adapter fitting to the most commonly used fiber optic cables, has an extremely high light carrier quality with good illumination in the complete operating field. The fiber optic channel is built to go inside the handle to increase the view during surgery avoiding that the fiber optic cable is disturbing the operation.



Blattgröße blade size	Artikelnummer order number
80 x 16 mm	260-316-081
80 x 24 mm	260-324-081
90 x 16 mm	260-316-091
90 x 24 mm	260-324-091
180 x 24 mm	260-324-181
190 x 16 mm	260-316-191
190 x 24 mm	260-324-191

mit Absaugung und Zähnen
with suction and teeth



Blattgröße blade size	Artikelnummer order number
80 x 16 mm	260-233-080
90 x 25 mm	260-247-090
110 x 16 mm	260-233-110
150 x 16 mm	260-234-150
150 x 25 mm	260-247-150
150 x 40 mm	260-248-150
180 x 16 mm	260-235-180
180 x 25 mm	260-247-180
220 x 16 mm	260-236-220
220 x 25 mm	260-247-220

mit Kaltlicht/with Fiber optic

SCULPO RETRAKTOREN
SCULPO RETRACTORS



mit Zähnen/with teeth


Blattgröße blade size	Artikelnummer order number
80 x 16 mm	260-080-160
90 x 20 mm	260-120-090
90 x 22 mm	260-122-090
90 x 25 mm	260-125-090
90 x 30 mm	260-137-090
120 x 22 mm	260-120-250
120 x 30 mm	260-120-300
125 x 25 mm	260-230-125
130 x 16 mm	260-034-130
130 x 16 mm	260-137-130
130 x 30 mm	260-130-130
150 x 36 mm	260-137-150
160 x 36 mm	260-136-162
180 x 25 mm	260-130-010
210 x 27 mm	260-130-210
220 x 36 mm	260-137-220

BRUSTRETRAKTOREN
BREASTRETRACTORS




11C

SCULPO RETRAKTOREN MIT KALTLICHT
SCULPO RETRACTORS WITH FIBER OPTIC




Blattgröße blade size	Blattwinkel blade angulation	Artikelnummer order number
80 x 25 mm		260-325-080
90 x 24 mm		260-134-090
90 x 30 mm		260-134-095
100 x 25 mm		260-135-100
100 x 25 mm	70°	260-135-101
125 x 25 mm		260-120-125
120 x 30 mm		260-130-120
130 x 16 mm		260-134-130
140 x 24 mm		260-213-141
150 x 25 mm		260-325-150
150 x 30 mm		260-130-150
150 x 36 mm		260-134-150
160 x 36 mm		260-136-161
190 x 36 mm		260-134-190



Blattgröße blade size	Blattwinkel blade angulation	Artikelnummer order number
70 x 16 mm		260-230-015
90 x 18 mm		260-230-018
90 x 22 mm		260-230-112
90 x 27 mm		260-230-017
100 x 25 mm		260-230-116
100 x 25 mm	70°	260-230-118
130 x 16 mm		260-230-019
150 x 16 mm		260-150-016
150 x 25 mm	70°	260-230-114
150 x 25 mm		260-230-014
150 x 30 mm		260-230-020
180 x 25 mm		260-230-016
220 x 27 mm		260-230-010

SCULPO RETRAKTOREN MIT ABSAUGUNG
SCULPO RETRACTORS WITH SUCTION



Blattgröße blade size	Besonderheiten specials	Artikelnummer order number
70 x 24 mm	Leichtmetall Handgriff/light metal handle	260-224-071
80 x 16 mm	Leichtmetall Handgriff/light metal handle	260-216-081
80 x 24 mm	Leichtmetall Handgriff/light metal handle	260-224-081
90 x 16 mm	Leichtmetall Handgriff/light metal handle	260-216-091
90 x 24 mm	Leichtmetall Handgriff/light metal handle	260-224-091
90 x 30 mm		260-119-300
150 x 36 mm		260-139-150
180 x 16 mm	Leichtmetall Handgriff/light metal handle	260-216-181
180 x 24 mm	Leichtmetall Handgriff/light metal handle	260-224-181
190 x 16 mm	Leichtmetall Handgriff/light metal handle	260-216-191

SCULPO RETRAKTOREN MIT ABSAUGUNG
SCULPO RETRACTORS WITH SUCTION




Blattgröße blade size	Blattwinkel blade angulation	Artikelnummer order number
80 x 16 mm		260-133-070
90 x 27 mm		260-129-270
120 x 14 mm		260-114-120
150 x 25 mm		260-150-250
180 x 16 mm	120° + 130°	260-130-216
180 x 25 mm	120° + 130°	260-130-225
180 x 25 mm		260-130-125
180 x 27 mm	120° + 130°	260-130-227
180 x 30 mm	120° + 150°	260-130-200

SCULPO RETRAKTOREN MIT KALTLICHT UND ABSAUGUNG
SCULPO RETRACTORS WITH FIBER OPTIC AND SUCTION




Die SCULPO Retraktoren überzeugen durch Ihre Anwendungsvielfalt und den griffigen Leichtmetall-Handgriff durch welchen wir unsere Kaltlichtkanäle auch in Kombination mit Absaugrohren legen können. In Zusammenarbeit mit vielen plastischen Chirurgen aus aller Welt setzt die RZ Medizintechnik GmbH neue Maßstäbe in der Mammachirurgie.

The SCULPO Retractors convince by the possibility of using them for a large variety of applications and by their very handy light metal handles, which are also supplied with a combination of fiber optic as well as suction. In collaboration with many plastic surgeons all over the world, RZ Medizintechnik sets new standards in mamma surgery.




Blattgröße blade size	Blattwinkel blade angulation	Artikelnummer order number
70 x 16 mm		260-216-071
70 x 25 mm		260-225-070
80 x 16 mm		260-216-082
80 x 24 mm		260-224-082
80 x 25 mm		260-135-008
90 x 15 mm		260-109-150
90 x 16 mm		260-216-092
90 x 22 mm		260-109-220
90 x 24 mm		260-109-240
90 x 24 mm		260-224-090
90 x 30 mm		260-109-300
100 x 24 mm	70°	260-224-191
100 x 25 mm		260-224-182
110 x 30 mm		260-109-110
120 x 36 mm		260-136-120
130 x 16 mm		260-216-132
150 x 16 mm		260-135-016
150 x 25 mm		260-125-150
150 x 30 mm		260-135-150
150 x 36 mm		260-136-150
160 x 36 mm		260-136-160
180 x 16 mm		260-180-016
180 x 25 mm		260-135-180
190 x 16 mm		260-216-192
190 x 24 mm		260-224-192
190 x 25 mm		260-135-190
190 x 30 mm		260-230-190
190 x 36 mm		260-136-190
190 x 50 mm		260-150-190
220 x 36 mm		260-136-220

SCULPO RETRAKTOREN MIT KALTLICHT UND ABSAUGUNG
SCULPO RETRACTORS WITH FIBER OPTIC AND SUCTION



Blattgröße blade size	Blattwinkel blade angulation	Besonderheiten specials	Artikelnummer order number
70 x 15 mm			260-107-150
70 x 20 mm			260-107-200
70 x 25 mm		mit isolierter Spitze/with insulated tip	260-370-025
90 x 16 mm			260-160-090
90 x 25 mm			260-190-240
90 x 27 mm			260-230-012
90 x 30 mm			260-190-300
100 x 25 mm			260-123-100
100 x 25 mm	70°		260-123-101
110 x 25 mm			260-125-110
120 x 25 mm			260-125-120
120 x 25 mm			260-244-020
150 x 16 mm		komplett beschichtetes Blatt/fully insulated blade	260-215-016
150 x 23 mm			260-123-150
150 x 25 mm			260-244-150
150 x 25 mm			260-142-150
150 x 25 mm	70°		260-134-151
150 x 25 mm	20°		260-135-151
150 x 30 mm			260-110-300
150 x 30 mm			260-230-150
150 x 36 mm			260-150-360
150 x 36 mm			260-243-150
180 x 16 mm			260-216-182
180 x 25 mm			260-244-180
180 x 25 mm			260-245-180
180 x 25 mm			260-136-001
180 x 30 mm			260-230-013
190 x 25 mm			260-230-011
190 x 36 mm			260-243-190
190 x 36 mm			260-136-019
220 x 25 mm			260-246-220
220 x 27 mm			260-136-000

SCULPO OPTIK RETRAKTOREN
SCULPO OPTIC RETRACTORS



Blattgröße blade size	Blattwinkel blade angulation	Besonderheiten specials	Artikelnummer order number
180 x 16 mm		für 4 mm Optik / for 4 mm optic	260-130-160
180 x 25 mm		für 4 mm Optik / for 4 mm optic	260-130-000
180 x 25 mm		für 10 mm Optik / for 10 mm optic	260-130-015
220 x 27 mm		für 4 mm Optik / for 4 mm optic	260-130-100
190 x 22 mm	35°	für 4 mm Optik / for 4 mm optic	260-130-194
190 x 22 mm	35°	für 4 mm Optik / for 4 mm optic	260-130-040
180 x 25 mm		für 5 mm Optik / for 5 mm optic	260-130-014