

Insulation and Continuity Tester

Have you checked your instrument insulation lately?

Our new 9-volt battery insulation and continuity tester is the perfect solution for locating potentially dangerous tears and cracks along the tube and handle of your Laparoscopic instruments.



#ATI-014


complete including accessories
and carrying case



#ATI-014 FC

complete including accessories,
carrying case and foot control
pedal

Insulation verification of




Laparoscopic electrodes, scissors
and forceps



Bayonet forceps

Continuity and short circuit test



Bipolar cables and
Laparoscopic bipolar
forceps

**2nd
Generation**

**Completely
Reusable**

For more information:

Jac-Cell Medic
105 Guthrie Avenue
Dorval (Québec)
Canada H9P 2P1

E-mail : mark@jaccell.com
Tel : (800) 720-0931
Fax : (514) 631-1351
Tel : (800) 667-6188 B.C. & Alberta