



Have you checked your instrument insulation lately?

The Original Portable Battery operated Insulation + Continuity Tester

Introducing The 4th Generation Insulation Tester with Footcontrol

New for 2016

Your Portable Insulation
Continuity Testing Device



#ATI-014AFC

For verifying:

- Insulated laparoscopic handles and tubes
- Insulated monopolar and bipolar bayonets
- Insulated leep vaginal speculums and forceps


For more information:

www.jacell.com

mark@jacell.com

New for 2015

3rd Generation
Microcontroller Driven
Completely Reusable Accessories



#ATI-021FC

We TEST at 2800 volts
DC (continuous tension)

Don't be fooled by
lower voltage testing
devices.



(800) 720-0931